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LENTI OFTALMICHE

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INSTRUCTIONS FOR USE, OPHTHALMIC LENSES

All lenses manufactured by DIVEL ITALIA S.p.A meet the requirements of EU Regulation 2017/745 about medical devices.

Optical lens performances follow the European standard EN ISO 14889 and the EN ISO 8980 series.

The Biological compatibility, the material inflammability and the mechanical strength have been verified in accordance with the EN ISO 14889. The transmittance specifications are verified following the test methods described in the ISO 8980-3 standard.

INSTRUCTIONS FOR USE

To ensure a longer life of the medical device we recommend:

- Be careful when you handle your glasses, don't put lenses in direct contact with hard surfaces;
- Wash your lenses only with water and neutral soap, without using aggressive products such as solvents, alcohol, acetone, etc., and using soft cloths to dry them;
- Don't clean the anti-reflective coating in dry conditions, especially when lenses are dusty and/or dirty;
- Don't expose your lenses to high temperatures, especially in summer. For example, don't leave your glasses onto the car dashboard.

Warnings

Uncut Finished Ophthalmic Lenses are provided to Wholesalers and Optician's laboratories. Although the easy edging process, it must be done only by specialized workers. The medical device (lens) is not designed for extreme conditions of use which could cause the damage of the device and could put the end user in danger.

If the lens must be tint, the technician in charge of the process must define the filter category and put it in declaration of conformity in compliance with EU Regulation 2017/745.

Recommendations

It is imperative to provide the original packaging and the readable label to the end-user. Neither the manufacturer nor the distributor shall be responsible in the event of the nonobservance of any of the instructions stated above.

In the case of critical incident involving the medical device produced by Divel Italia, it must be reported to the manufacturer and to the relevant authority of the member state which it belongs.

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