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INSTRUCTIONS FOR USE,

OPHTHALMIC LENSES

Ophthalmic lenses are optical devices designed to correct visual defects such as myopia, hyperopia, astigmatism, or presbyopia.

All lenses manufactured by **DIVEL ITALIA S.p.A** meet the requirements of EU Regulation 2017/745 on medical devices.

Optical lens performances follow the European standard EN ISO 14889 and the EN ISO 8980 series.

Physiological compatibility, flammability of the materials used, and the mechanical strength of the lenses have been verified according to EN ISO 14889. Absorption values are determined using the testing methods described in ISO 8980-3.

WARNINGS

Ophthalmic lenses are intended for wholesalers and optical laboratories. The lens edging and mounting process, although simple, must be done exclusively by qualified personnel. Improper use or processing with unsuitable tools may compromise optical quality.

The medical device (lens) is not designed for extreme conditions of use which could cause the damage of the device and could put the end user in danger.

If the lens must be tint, the technician in charge of the process must define the filter category and put it in declaration of conformity in compliance with EU Regulation 2017/745.

INSTRUCTIONS FOR USE

- Always verify that optical specifications (power, axis, prism) correspond to the prescription before processing the lenses.
- Wear microfiber gloves to prevent scratches or fingerprints.
- Use only calibrated and clean cutting and polishing tools.
- Avoid overheating during edging to prevent surface deformation and/or damage to the anti-reflective coating.
- Ensure that the edges are smooth and free of micro-fractures.

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TEL.
+39 051 721651

FAX
+39 051 721860

- Clean lenses only after removing any processing residues.
- Before mounting lenses into the frame, checking horizontal and vertical centration.
- Before delivering the glasses to the end user, thorough verification must be carried out to ensure that they conform to the medical prescription and the user's individual optical parameters.

CLEANING AND MAINTENANCE

To ensure a longer life of the medical it is recommended to

- Be careful when you handle your glasses, don't put lenses in direct contact with hard surfaces;
- Wash your lenses only with water and neutral soap, without using aggressive products such as solvents, alcohol, acetone, etc., and using soft cloths to dry them;
- Don't clean the anti-reflective coating in dry conditions, especially when lenses are dusty and/or dirty;
- Don't expose your lenses to high temperatures, especially in summer. For example, don't leave your glasses onto the car dashboard.

RECOMMENDATIONS

It is imperative to provide the original packaging and the readable label to the end-user. Neither the manufacturer nor the distributor shall be responsible in the event of the nonobservance of any of the instructions stated above.

In the case of critical incident involving the medical device produced by Divel Italia, it must be reported to the manufacturer and to the relevant authority of the member state which it belongs.

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